

>PRIMERS (all primers are listed 5'---->3')

>

>Exon 11: anti-sense

>ttg gtt ctc ttg aac tct gc

>

>Neo 2: anti-sense

>act gca tct gcg tgt tcg aa

>

>Intron: sense; common primer

>taa gag taa tct tcc aga gc

>

>REACTION_____

>

>10.0 uL : DNA (20 ng) of genomic per reaction (I use ~0.3 uL of DNA)

> 5.0 uL : 10xBuffer

> 0.4 uL : 25 mM dNTP

> 2.5 uL : [20 uM] Exon 11

> 2.5 uL : [20 uM] Neo 2

> 5.0 uL : [20 uM] Intron 5

> 0.5 uL : Taq (Fisher)

>26.5 uL : H2O

>-----

>50.0 uL Total

>

>CYCLING_____

>

>step 1: 94, 3 min

>step 2: 1.5 degrees/sec to 55

>step 3: 55, 45 sec

>step 4: 72, 1 min

>step 5: 94, 45 sec

>step 6: 1.5 degrees/sec to 55

>step 7: 55, 45 sec

>step 8: 72, 1 min

>step 9: goto step 5, 35 cycles

>step 10: 72, 5 min

>step 11: 4, for ever

>GEL_____

>

>run 25 uL on a 1.0% agarose gel for ~1:30 hours at 90V.

>

>bands: (exon 11, neo 2, intron 5)

>wild-type: 426 bp

>heterozyg: 426 bp + 600 bp

>knock-out: 600 bp